

Status Report: Peroxides and Formaldehyde R/V Atlantis, Spring 2010

Brian Heikes
Graduate School of Oceanography
University of Rhode Island
bheikes@gso.uri.edu



11-13 January 2011



Species Measured

- Hydrogen Peroxide, H_2O_2
- Methylhydroperoxide, CH_3OOH
- Formaldehyde, CH_2O aka HCHO
- Photochemical significance
 - VOCs → CH_2O → CO → CO_2
 - CH_2O → HOx → O_3
 - HOx → H_2O_2 → HOx
 - CH_4 and some OVOC → CH_3OO → CH_3OOH
 - CH_3OO → CH_3OOH → CH_2O and HOx

Climate Connections

CO₂ (+1.6), CH₄ (+0.5), O₃ (+0.4)
and sulfate (-0.4)

- VOCs and OVOCs → CH₂O → CO.
- CH₄ → CH₂O
 - and CH₄ → CH₃OOH.
- Peroxides and formaldehyde are directly relatable to ozone.
- Heterogeneous sulfur dioxide oxidation.

Climate Connections


adapted from Ricky Rood (UMich),
“blog” Weather Underground, June 2010

Global radiative forcing
due to emissions of
aerosols and precursors
changes from
1750-2005

Black Carbon - - - -
Positive Radiative
Forcing

Sulfate and Organic
Carbon - - - - - Negative
Radiative Forcing

Cloud Albedo Effect
-- Negative Radiative
Forcing



Formaldehyde Method - I

- Gas-to-aqueous collection
- Enzyme catalyzed derivatization
- Fluorescence product detection

- ~80% collection efficiency in CalNex
- Gas and aqueous stds.
- LN2 used for carrier gas and blanks

Formaldehyde Status

- Temporal resolution 2 minutes.
- Significant lag, ~ 25 minutes.
- LOD, 0.05 to 0.10 ppb.
- C-tainer temp effects and peroxide “cross-talk”.
 - Need to remove data under conditions of “rapid” change.
 - Need to remove data from periods of peroxide repair.
- Near final data reduction; ~ 1 month to go.


Peroxide Method - I

- Gas-to-aqueous collection.
 - HPLC separation of hydroperoxides.
 - Enzyme catalyzed derivatization.
 - Fluorescence detection.
-
- 100 / 70% collection of H_2O_2 / CH_3OOH .
 - Gas and aqueous stds.
 - LN_2 used for carrier gas and blanks.


Peroxide Method - Status

- Nearly every component failed at some point in the deployment.
- 30-sec sample every 4 minutes.
- Precision
 - H₂O₂, 0.01 to 0.05 ppb
 - CH₃OOH 0.03 to 0.1 ppb
- Hand analysis of chromatograms - time consuming.
 - Four more months to an archival product.

Example Data: 5/18-19/2010



Example Data June 1-2, 2010



Monterey to San Francisco Bays

Summary

- Formaldehyde data in ~1 month.
- Peroxide data in ~4 months.
- Dynamic photochemical interplay between peroxides and between peroxides and formaldehyde.
- Will examine BL height, cloud and ocean surface and hydrogen peroxide.

It's a small world at times;
sometimes even smaller.

